

MULTISAMPLE 120 PHANTOM



The MultiSample120 phantom is a 120mm diameter cylindrical phantom, capable of holding up to eight industry standard Falcon BD 15ml polypropylene tubes. The phantom can be supplied empty for filling by the user, or can be filled with solutions to the customer's specification.

Uses include:

- Quantitative T1/T2 relaxometry.
- Metabolic imaging and spectroscopy.
- Multi-nuclear imaging and spectroscopy.

SPECIFICATIONS

Diameter	120mm (internal volume) 130mm (outer cylinder) 155mm (widest point)*
Length	165mm*
Overall Dimensions (L/W/H)	155mm x 155mm x 165mm*
Materials	Acrylic (body, lid and) Nylon (fasteners) Nitrile (gasket seal) HDPE (plug)
Sample containers	8x Falcon BD 15ml Polypropylene Tubes Easy snap-in design, tubes can be easily changed without the need for tools.
Contents	Can be supplied empty or filled to customer specification
Packaging	Foam padded box for storage and transport
Optional Extras	Connection fittings for a heated water system

*Dimensions will be adjusted to ensure fit with the specified MRI coils.

Specifications/images are provided as a guide only and may be subject to change.

Gold Standard Phantoms Limited

Registered address Unit 103 Belgravia Workshops, 159-163 Marlborough Road, London, N19 4NF
Registered in England and Wales. Company number: 9342804. VAT REG: GB211807049